

The Welch Allyn Green Series™ 600 Minor Procedure Light



Disposable Sheath
Available—Use #52630



Easy Power
Connector



3 LEDs



Touchless On/Off
with Dimming Control

Easy Assembly



WelchAllyn®

Advancing Frontline Care™

Green Series™ Medical Exam Lights



Whiter, Brighter, and Greener

The GS 600 Minor Procedure Light features:

- Three LEDs and a wider head, designed specifically for demanding viewing conditions
- Uniform spot and superior brightness for minor procedures and exams in any healthcare environment
- Less energy consumption
- No bulb replacement



IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS

- PRIMARY CARE CLINICS
- OB/GYN
- OUTPATIENT FACILITIES
- DERMATOLOGY
- HOSPITAL—ER AND ICU

MOUNTING OPTIONS AVAILABLE:



Table/Wall Mount

Mobile Stand

Light Specifications	# of LEDs	3 White LEDs
	Colour Temperature	5,500K
	Light Intensity @ Working Distance	50,000 Lux @ 30cm (12in) 45,000Lux @ 50cm (20in)
	Spot Size @ Working Distance	18cm @ 30cm (12in) 15cm @ 50cm (20in)
Physical Dimensions	Reach Mobile Stand	79cm (30,1in)
	Reach Wall	79cm (30,1in)
	Reach Ceiling	n/a
	Ceiling up/down	n/a
	Luminaire (Height/Width/Depth)	34,5 / 19,2 / 8,6 cm (13,58 / 7,56 / 3,39in)
	Control Box (Height/Width/Depth)	27,5 / 10,0 / 8,0 cm (10,83 / 3,94 / 3,15in)
	Mobile Stand (Height / Diameter)	61cm / 53cm (24,02in / 20,87in)
	Total Weight Table/Wall mount	2,93 kg (6,45 lbs)
	Total Weight Mobile Stand	7,87 kg (17,35 lbs)
Electrical	Power Supply	Input: 100-240 V~ 50Hz-60Hz, 0,65A max
	Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2,5 m) long
	Control Functions	Touch less On/Off with Dimming
Misc Data	Mounting Options	Mobile Stand, Table/Wall, Rail
	Kind of Extension	Gooseneck
	Accessory Options	Disposable sheath: P/N 52630 Length: 22,86 cm (9in)
	Lifetime	50,000 hours
	Warranty	5-year limited warranty

To place an order, contact your appropriate Customer Service Department.



Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220, U.S.A.
Telephone: +1.315.685.4602
+1.800.535.6663
Fax: +1.315.685.0054

Welch Allyn European Customer Service: +353 46 90 67790
Australia: +61 29 638 3000
China: +86 21 6327 9631
France: +33 160 09 33 66
France (Customer Service): +33 155 69 58 49
Germany: +49 6950 985 132
Spain: +34 912 990 376
Japan: +81 3 3219 0071
India: +91 44 4353 5724

Latin America (Miami office): +1 305 669 9003
Malaysia: +603 78753341
Netherlands: +31 202 061 360
Singapore: +65 6419 8100
Sweden: +46 042 33 02 20
United Kingdom: +44 207 365 6780

www.welchallyn.com/greenseris

WelchAllyn®

Advancing Frontline Care™

