

The Welch Allyn Green Series™ Exam Light IV



Highly Flexible Light Pipe—
Twice as long as conventional pipes

1 LED—3X the typical light
output; pinpoint accuracy;
true tissue color rendition



**Disposable Sheath
Available—**Use existing
disposable sheath (#52640);
also for use with Exam Light
II and III



**Easy Power
Connector**



Touchless On/Off



Easy Assembly

WelchAllyn®

Advancing Frontline Care™

Green Series™ Medical Exam Lights



Whiter, Brighter, and Greener

The GS Exam Light IV features:

- A compact, unobtrusive head, ensuring pinpoint accuracy
- A whiter, brighter light than halogen, so you see true tissue color rendition during an exam
- Less energy consumption
- No bulb replacement



Light Specifications	# of LEDs	1 White LED
	Colour Temperature	5,500K
	Light Intensity @ Working Distance	20,000 Lux @ 40cm (16in)
	Spot Size @ Working Distance	10-25cm (4-10in) @ 40cm (16in)
Physical Dimensions	Reach Mobile Stand	126cm (49,61in)
	Reach Wall	126cm (49,61in)
	Reach Ceiling	n/a
	Ceiling up/down	n/a
	Luminaire (Height/Width/Depth)	Length: 13,1cm (5,16in) Diam.: 4cm (1,58in)
	Control Box (Height/Width/Depth)	27,5 / 10,0 / 8,0 cm (10,83 / 3,94 / 3,15in)
	Mobile Stand (Height / Diameter)	61cm / 53cm (24,02in / 20,87in)
	Total Weight Table/Wall mount	2,4 kg (5,30 lbs)
	Total Weight Mobile Stand	5,63 kg (12,40 lbs)
Electrical	Power Supply	Input: 100-240 V~ 50Hz-60Hz, 0,30A max
	Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2,5 m) long
	Control Functions	Touch less On/Off
Misc Data	Mounting Options	Mobile Stand, Table/Wall, Rail
	Kind of Extension	Gooseneck
	Accessory Options	Disposable sheath: P/N 52640 Length: 129.54 cm (51in)
	Lifetime	50,000 hours
	Warranty	5-year limited warranty

IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- OB/GYN
- FAMILY PRACTICE
- PLASTIC SURGERY
- DERMATOLOGY
- EMERGENCY ROOMS
- OUTPATIENT SURGERY
- LABOR AND DELIVERY SETTINGS
- SURGI-CENTERS

MOUNTING OPTIONS AVAILABLE:



Table Mount



Wall Mount



Mobile Stand

To place an order, contact your appropriate Customer Service Department.



www.welchallyn.com/greenserier

Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220, U.S.A.
Telephone: +1.315.685.4602
+1.800.535.6663
Fax: +1.315.685.0054

Welch Allyn European Customer Service: +353 46 90 67790
Australia: +61 29 638 3000
China: +86 21 6327 9631
France: +33 160 09 33 66
France (Customer Service): +33 155 69 58 49
Germany: +49 6950 985 132
Spain: +34 912 990 376
Japan: +81 3 3219 0071
India: +91 44 4353 5724

Latin America (Miami office): +1 305 669 9003
Malaysia: +603 78753341
Netherlands: +31 202 061 360
Singapore: +65 6419 8100
Sweden: +46 042 33 02 20
United Kingdom: +44 207 365 6780

WelchAllyn®

Advancing Frontline Care™

